

Single Patient Use Tourniquet

The Gold Standard of Single Patient Use Tourniquets

Single Patient Use

Designed for effective infection prevention

Single-handed engagement and release

Safe and effective in management of blood sampling and IV cannulation

Quick release and pull through fastening

Remains secure and in position during blood drawing

Super soft material for increased patient comfort


Unique formulated material designed to be strong but soft next to the skin

Integral dispenser pack

Re-sealable packaging for increased hygiene

For Paediatric & Adult Use


One size fits ages 2 months and above



Single Patient Use Tourniquet

The tourniSTRETCH® from ASep Healthcare was developed by Healthcare Professionals to provide users with a premium Single Patient Use tourniquet for everyday use, challenging blood sampling and venepuncture. Made from a customised material, it provides users with strong, secure stretch for optimal tensioning whilst maintaining a soft flexible feel for improved patient experience. It's intuitive pull through system, easy re-adjustment and quick release provides professionals with easy application and assists in single-handed operation when blood sampling.

Instructions



Pass TourniSTRETCH round the limb peg side up

Feed tail through buckle and tension accordingly

Pull tail back over the peg to secure in position

Product Code: 4345ITPEO/ASEP25



Integral Resealable dispenser pack

Contains 25 Single Patient Use tourniquets

Single pack option for clean storage for each patient



TourniSTRETCH is a Trade Mark of ASep Healthcare Ltd

TourniSTRETCH is distributed by ASep Healthcare Ltd

rev 2.0 2020-01

T: +44 (0)151 319 2963

support@asephealthcare.com

18 Tapton Way, Wavertree Business Village,
Liverpool L13 1DA

www.asephealthcare.com